

Isotemp Water Baths



Consistent Temperature Control
for Reliable Results



37.0°C

Isotemp[®] GPD 15D



37.0°C

Isotemp



Fisher Scientific Isotemp General Purpose Water Baths

Ideal for industrial, clinical pharmaceutical or biomedical applications, Fisher Scientific™ Isotemp™ general purpose water baths provide powerful technology in a compact design. Available in a variety of sizes to meet a wide range of application requirements, Isotemp water baths provide consistent temperature for reliable results in an easy-to-use operation.

Features:

- » New timer for worry free operation
- » New addition of three presets for quick set up
- » Advanced temperature controller with LCD display shows set-point and actual temperature
- » Global voltage – each bath works anywhere in the world exactly the same
- » Corrosion-resistant exterior with easy-to-clean enamel coating
- » Display shows °C or °F and is readable from a distance
- » Open tank design makes the bath easy to clean
- » Adjustable over temperature protection
- » Low fluid start protection capability
- » CE and UL listed
- » Audible alarms for timer, low fluid detection, over temperature and safety
- » Rounded reservoir corner and sample tray with finger holes allow for easy removal and cleaning
- » Hinged polycarbonate cover that lifts to a 90° stay-open position or can be removed completely
- » Baths come with diffuser tray and gable cover



Fisher Scientific Isotemp General Purpose Water Bath Specifications

Model	Cat. No.	Chamber Capacity	Temperature Range	Temperature Stability/Uniformity @ 70°C	Work Area (L x W x H)	Overall Dimensions Without Cover (L x W x H)	Global Voltage	Heater Output [†]
GPD 02	15355877	2L	Amb. to 90°C	±0.1°C / ±0.2°C	5.4 x 6.1 x 5.9" (138 x 155 x 150mm)	9.1 x 7.8 x 9.2" (230 x 199 x 233mm)	100-115V/200-230V, 50/60Hz	200W
GPD 02S	15365877	2L (Shallow)	Amb. to 100°C	±0.1°C / ±0.2°C	6 x 11.8 x 2.6" (153 x 300 x 65mm)	9.7 x 14 x 9.1" (246 x 355 x 232mm)	100-115V/200-230V, 50/60Hz	300W
GPD 05	15375877	5L	Amb. to 100°C	±0.1°C / ±0.2°C	6.1 x 11.8 x 5.9" (154 x 300 x 150mm)	9.7 x 14 x 9.1" (246 x 355 x 232mm)	100-115V/200-230V, 50/60Hz	300W
GPD 10	15385877	10L	Amb. to 100°C	±0.1°C / ±0.2°C	11.9 x 13 x 5.9" (301 x 330 x 150mm)	15.5 x 15.1 x 9.2" (393 x 383 x 233mm)	100-115V/200-230V, 50/60Hz	800W
GPD 20	15395877	20L	Amb. to 100°C	±0.1°C / ±0.2°C	11.7 x 19.7 x 5.9" (297 x 500 x 150mm)	15.4 x 21.8 x 9.2" (392 x 555 x 233mm)	100-115V/200-230V, 50/60Hz	1200W
GPD 28	15305887	28L	Amb. to 100°C	±0.1°C / ±0.2°C	11.7 x 19.7 x 7.9" (297 x 500 x 200mm)	15.4 x 21.8 x 11.1" (392 x 555 x 282mm)	100-115V/200-230V, 50/60Hz	1200W
GPD 15D	15315887	5L & 10L (Dual)	Amb. to 100°C	±0.1°C / ±0.2°C	See FSGPD05 & FSGPD10	15.4 x 23.1 x 9.2" (392 x 587 x 233mm)	100-115V/200-230V, 50/60Hz	300W & 800W

[†]Heater Output at 120V and 240V

Accessories for Isotemp General Purpose Water Baths

	15355877	15365877	15375877	15385877	15395877	15305887	15315887
Stainless Steel Cover	15385887	15305897	15305897	15315897	15325897	15325897	15305897 & 15315897
Concentric Ring Cover	-	-	-	15355897	15365897	-	-
Replacement Acrylic Cover	15345887	15355887	15355887	15365887	15375887	15375887	15355887 & 15365887
Stainless Steel Petri Dish Rack	-	-	-	-	15385897	15385897	-
Stainless Steel Test Tube Rack	-	-	-	-	15395897	15395897	-
Hand Pump	15365927	15365927	15365927	15365927	15365927	15365927	15365927



Fisher Scientific Isotemp Shaking Water Baths



With a capacity up to 27 liters (7 gallons) Fisher Scientific™ Isotemp™ shaking water baths have intuitive digital control and are easy-to-use. These versatile shaking water baths can be used across a wide array of applications from life science and QA/QC. Count on the latest in technology, features and options to deliver the ideal combination of performance, ease of use and reliability.

Features:

- » New hinged stainless steel cover that lifts to a 90°C stay-open position or can be removed completely
- » Accurate temperature control and quiet reciprocal shaking action
- » Adjustable shaking speed from 30 to 200 oscillations/min (OPM)
- » Timer for worry free operation
- » Three presets for quick set up of temperature and shaking speed
- » Advanced temperature controller with LCD display shows set-point and actual temperature
- » Global voltage – each bath works anywhere in the world exactly the same
- » Corrosion-resistant exterior with easy-to-clean enamel coating
- » Display shows °C or °F is readable from a distance
- » Open tank design makes the bath easy to clean
- » Adjustable over temperature protection
- » Low fluid start protection capability
- » CE and UL Listed
- » Audible alarms for timer, low fluid detection, over temperature and safety
- » Baths come with shaking tray and gable cover



Fisher Scientific Isotemp Shaking Water Bath Specifications

Model	Cat. No.	Chamber Capacity	Temperature Range	Temperature Stability/Uniformity @ 70°C	Work Area (Removable Tray) (L x W x H)*	Overall Dimensions Without Cover (L x W x H)	Global Voltage	Heater Output†
SWB 15	15325887	15 Liter	Amb. + 5°C to 100°C	±0.1°C / ±0.05°C	11.5 x 12 x 6.5" (292 x 305 x 165mm)	15.5 x 24.9 x 9.8" (394 x 632 x 249mm)	100-115V/200-230V, 50/60Hz	1200W
SWB 27	15335887	27 Liter	Amb. + 5°C to 100°C	±0.1°C / ±0.05°C	11.5 x 24 x 6.5" (292 x 610 x 165mm)	15.5 x 36.9 x 9.8" (394 x 938 x 249mm)	100-115V/200-230V, 50/60Hz	1500W

*Fits in Bath Opening (L x W x H). On 15L Shaking Baths: 12 x 15.3 x 7.6" (305 x 387 x 193mm). On 27L Shaking Bath: 12 x 27.3 x 7.6" (305 x 692 x 193mm).

†Heater Output at 120V and 240V.



Accessories for Isotemp Shaking Water Baths

	15325887	15335887		15325887	15335887		15325887	15335887
0.5mL Microfuge Tube Rack*	15305917	15305917	Test Tube Clip 13mm-25mm*	15355917	15355917	Flask Clips – 250mL*	15395917	15395917
1.0mL Microfuge Tube Rack*	15315917	15315917	Fit in bath	20	48	Fit in bath	6	14
Water Regulation Kit	15355927	15355927	Flask Clips – 25mL*	15365917	15365917	Flask Clips – 500mL*	15305927	15305927
Test Tube Tray 13-25mm	15325917	15325917	Fit in bath	20	48	Fit in bath	4	12
Flask Tray 25mL	15335917	15335917	Flask Clips – 50mL*	15375917	15375917	Flask Clips – 1000mL*	15315927	15315927
Flask Tray 50mL	15345917	15345917	Fit in bath	15	36	Fit in bath	2	5
High Wall Tray – Small	15325927	-	Flask Clips – 125mL*	15385917	15385917	Fasteners (packs of 25)	15345927	15345927
High Wall Tray – Large	-	15335927	Fit in bath	9	24			

*Need Fasteners

Isotemp Water Baths

Consistent Temperature Control
for Reliable Results



© 2015 Thermo Fisher Scientific Inc. All rights reserved.
Trademarks used are owned as indicated at www.fishersci.com/trademarks.

Austria: (0)800-20 88 40 Belgium: +32 (0)56 260 260 Denmark: +45 70 27 99 20
Germany: +49 (0)2304 9325 Ireland: +353 (0)1 885 5854 Italy: +39 02 950 59 478
Finland: +358 (0)9 8027 6280 France: +33 (0)3 88 67 14 14 Netherlands: +31 (0)20 487 70 00
Norway: +47 22 95 59 59 Portugal: +351 21 425 33 50 Spain: +34 902 239 303

 **Fisher
Scientific**

A Thermo Fisher Scientific Brand

EU_EN_2015_11475